



## Tourniquet Systems

Tourniquet devices with cuffs and accessories

# / Tourniquet Touch TT20 / TT15



The Tourniquet Touch is an electrically operated surgical tourniquet. It regulates the pressure of a tourniquet cuff which temporarily occludes the blood flow of a patient's upper or lower extremity in order to obtain a bloodless field. User friendly innovations such as fast choice buttons and the specific alarm system simplify the use of the device and ensure patient safety. The system interface is reliable and intuitive and allows easy maintenance for the biomedical engineer (calibration, self test, leak test). The integrated lithium-ion battery allows a backup runtime of approx. 8 h.

## **Tourniquet Touch TT20 with Two Cuff Channels**

For use with single cuff, two single cuffs or double cuff.

## **Tourniquet Touch TT15 with Cuff Channel and Irrigation Channel**

For use with single cuff and pressure infusion cuff.

Order  
information

**Tourniquet Touch TT20** / With 3.0 m Coil Connecting Tubing (blue and red), 100 - 240 VAC, height 186 mm, width 263 mm, depth 226 mm, 4.5 kg

REF	Box
-----	-----

01-20-000	1
-----------	---

**Tourniquet Touch TT15** / With 3.0 m Coil Connecting Tubing (blue and black), 100 - 240 VAC, height 186 mm, width 263 mm, depth 226 mm, 4.5 kg

REF	Box
-----	-----

01-15-000	1
-----------	---

**Mobile Stand with Basket** / 4 castors with castor lock (2 include ESD), height 939 mm, width 400 mm, depth 400 mm, 8.2 kg

REF	Box
-----	-----

01-00-100	1
-----------	---

Additional  
information

Application video



Brochure

# / Tourniquet Dispo Cuff

Tourniquet cuff, for single use

Tourniquet Dispo Cuff is a sterile, single use tourniquet cuff especially suitable for contaminated procedures, preventing against HAI (Hospital-acquired infection). The cuff is available in straight shape (for small or midsize upper extremity) and contour shape (mainly for lower conical shaped extremity) and is provided in different sizes, which are colour coded for optimal cuff selection.



Order  
information**Tourniquet Dispo Cuff / Single cuff, for single use, sterile**

Length	Width	Colour	Single Port	Dual Port	Box
20 cm	5 cm	-	REF 20-34-700SLZ-1	REF 20-36-700SLZ-1	10
30 cm	6 cm	■	REF 20-34-710SLZ-1	REF 20-36-710SLZ-1	10
35 cm	8 cm	■	REF 20-34-711SLZ-1	REF 20-36-711SLZ-1	10
46 cm	10 cm	■	REF 20-34-712SLZ-1	REF 20-36-712SLZ-1	10
46 cm	12 cm	■	REF 20-34-715SLZ-1*	REF 20-36-715SLZ-1*	10
61 cm	12 cm	■	REF 20-34-722SLZ-1*	REF 20-36-722SLZ-1*	10
76 cm	12 cm	■	REF 20-34-727SLZ-1*	REF 20-36-727SLZ-1*	10
86 cm	12 cm	■	REF 20-34-728SLZ-1*	REF 20-36-728SLZ-1*	10
107 cm	12 cm	■	REF 20-34-729SLZ-1*	REF 20-36-729SLZ-1*	10

**Tourniquet Dispo Cuff / Double cuff, for single use, sterile**

Length	Width	Colour	REF	Box
30 cm	12 cm	■	20-30-710SLZ-1	10
46 cm	15 cm	■	20-30-712SLZ-1	10
61 cm	15 cm	■	20-30-722SLZ-1	10

\*Contour shape

**Extension Tubing / For single use, sterile**

Length	REF	Box
100 cm	20-39-100-1	5

Additional  
information

Application arm



Application leg












Brochure






Order  
information

## Tourniquet Wipe Cuff / Single cuff, reusable

Length	Width	Colour	REF	Box
20 cm	7 cm		20-75-700	1
30 cm	8 cm		20-75-710	1
35 cm	8 cm		20-75-711	1
46 cm	10 cm		20-75-712	1
46 cm	12 cm		20-75-715*	1
61 cm	12 cm		20-75-722*	1
76 cm	12 cm		20-75-727*	1
86 cm	12 cm		20-75-728*	1
107 cm	12 cm		20-75-729*	1

## Tourniquet Wipe Cuff / Double cuff, reusable

Length	Width	Colour	REF	Box
30 cm	12 cm		20-77-710	1
46 cm	15 cm		20-77-712	1
61 cm	15 cm		20-77-722	1

\*Contour shape

Additional  
information

Application arm



Application leg




Brochure

# / Tourniquet Silicone Cuff

Tourniquet cuff, reusable through automatic processing










Tourniquet Silicone Cuff is a reusable tourniquet cuff made of durable, heavy-duty material and is suitable for repeated use after sterilization by steam autoclave at a temperature of 134°C for 5 minutes. The cuff is available in straight shape (for small or midsize upper extremity) and contour shape (mainly for lower conical shaped extremity) and is provided in different sizes, which are colour coded for optimal cuff selection. UDI code and serial number allow product traceability and tracking of product usage.








Order  
information**Tourniquet Silicone Cuff /** Single cuff, reusable, autoclavable

Length	Width	Colour	REF	Box
20 cm	2.5 cm	-	20-64-700	1
30 cm	4.5 cm		20-64-710	1
35 cm	7 cm		20-64-611	1
46 cm	8.5 cm		20-64-612	1
46 cm	11 cm		20-64-512*	1
61 cm	11 cm		20-64-522*	1
76 cm	11 cm		20-64-527*	1
86 cm	11 cm		20-64-528*	1

**Tourniquet Silicone Cuff /** Double cuff, reusable, autoclavable

Length	Width	Colour	REF	Box
30 cm	8 cm		20-60-710	1
46 cm	12.5 cm		20-60-712	1
61 cm	12.5 cm		20-60-722	1

\*Contour shape

**Sealing Cap /** For Tourniquet Silicone Cuff, reusable

REF	Box
22-60-111	5

Additional  
information

Application arm



Application leg



Brochure

# / Measuring Tape


For selection of VBM tourniquet cuffs

- Colour coded, for easy identification of cuff selection

Order  
information

**Measuring Tape** / For single use

REF	Box
20-95-222	50



# / Tourniquet Sleeve


For patient protection

- Protects soft tissues underneath the cuff
- Reduces wrinkling, pinching and shearing of soft tissues

Order  
information

**Tourniquet Sleeve** / For single use

Extremity	Colour	Clean Packed	Sterile Packed	Box
13 - 18 cm	<input type="checkbox"/>	REF 23-00-700	REF 23-00-700-1	10
18 - 27 cm	<input type="checkbox"/>	REF 23-00-710	REF 23-00-710-1	10
25 - 32 cm	<input type="checkbox"/>	REF 23-00-711	REF 23-00-711-1	10
32 - 43 cm	<input type="checkbox"/>	REF 23-00-712	REF 23-00-712-1	10
43 - 57 cm	<input type="checkbox"/>	REF 23-00-722	REF 23-00-722-1	10
57 - 69 cm	<input type="checkbox"/>	REF 23-00-727	REF 23-00-727-1	10
69 - 79 cm	<input type="checkbox"/>	REF 23-00-728	REF 23-00-728-1	10
79 - 100 cm	<input type="checkbox"/>	REF 23-00-729	REF 23-00-729-1	10



# / Tourniquet Drape



For patient protection

- Isolates the cuff from the sterile field
- Seals against fluid accumulation underneath the cuff to protect the patient from burns associated with skin preparation

Order  
information

**Tourniquet Drape** / For single use

Extremity	Clean Packed	Box	Sterile Packed	Box
18 - 45 cm	REF 24-00-100	40	REF 24-00-100-1	10
40 - 80 cm	REF 24-00-200	40	REF 24-00-200-1	10



# / Cuff Cover

For cuff protection

- Keeps the cuff free from soiling

Order  
information

**Cuff Cover** / For single use

Extremity	Size	Clean Packed	Sterile Packed	Box
18-43 cm	Small	REF 20-11-111	REF 20-11-111-1	50
40-61cm	Medium	REF 20-11-333	REF 20-11-333-1	50
58-81 cm	Large	REF 20-11-222	REF 20-11-222-1	50



# / Esmarch Bandages

For exsanguination

- Available in 4 sizes
- Made of synthetic rubber and powder-free (talcum-free) to avoid possible allergic reactions
- Textured surface for a good grip
- Elastic properties for easy wrapping
- Robust and tear resistant


Order  
information

**Esmarch Bandage /** For single use, sterile packed, colour: light blue

Width	Length	REF	Box
6 cm	365 cm	20-21-060-1	10
8 cm	365 cm	20-21-080-1	10
10 cm	365 cm	20-21-100-1	10
15 cm	365 cm	20-21-150-1	10

**Esmarch Bandage /** Reusable, autoclavable, colour: dark blue

Width	Length	REF	Box
6 cm	500 cm	20-20-060	5
8 cm	500 cm	20-20-080	5
10 cm	500 cm	20-20-100	5
15 cm	500 cm	20-20-150	5








# / Roll-On Cuff


For exsanguination

- Allow fast and easy exsanguination within seconds
- Inflatable exsanguination system
- 100% silicone medical grade
- Quick and efficient procedure
- Full range of sizes available

Order  
information

**Roll-On Cuff / Reusable**

Size	Extremity	Colour	REF	Box
2	26 - 33 cm		21-95-712	1
3	31 - 40 cm		21-95-719	1
4	38 - 49 cm		21-95-722	1
5	47 - 60 cm		21-95-727	1
6	58 - 70 cm		21-95-729	1



# / Pressure Infusion Cuffs

Administration of large fluid volumes for irrigation of the operation field



**Classic Infusor**

With zipper, allows insertion of different fluid bag configurations.



**Clear Infusor**

Clear 360° view on the fluid bag and fluid level.



**Dispo Infusor**

With hook, for easy fixation of most fluid bag configurations.



Order  
information**Classic Infusor** / Reusable

Version	3000 / 5000 ml	Box
With hand inflator, three-way stopcock and manometer	REF 58-08-300	1
With Luer Lock connection	REF 58-01-300	1



Application video



Brochure

500 ml and 1000 ml available on request

**Clear Infusor** / Reusable

Version	3000 ml	Box
With hand inflator, three-way stopcock and manometer	REF 57-08-300	1
With Luer Lock connection	REF 57-01-300	1



Application video



Brochure

500 ml and 1000 ml available on request

**Dispo Infusor** / For single patient use

Version	3000 ml	Box
With hand inflator, three-way stopcock and colour coded pressure indicator	REF 56-02-300	5
With Luer Lock connection	REF 56-01-300	5



Application video



Brochure

500 ml and 1000 ml available on request

# / Manual Tourniquet



The manual tourniquet regulates the pressure of a tourniquet cuff which temporarily occludes the blood flow of a patient's upper or lower extremity in order to create a bloodless field.

The hand inflator with manometer allows an autonomous use without electricity and without central air.

Order  
information

**Manual Tourniquet** / Hand inflator with manometer, pressure range max. 600 mmHg

With Blue Coil Connecting Tubing and Positive Locking Connector	With Red Coil Connecting Tubing and Positive Locking Connector	Box
REF 20-18-601	REF 20-18-602	1



Additional  
information



Find our full range of products on [www.vbm-medical.de](http://www.vbm-medical.de).

Here you will find our media centre with additional information such as application videos, brochures, flyer etc.

The medical devices in this advertising material are manufactured without the use of natural rubber latex, unless otherwise specified.

The medical devices in this advertising material do not contain phthalates which require labelling according to CLP Regulation (EC) 1272/2008.

**VBM** **Medizintechnik GmbH**

Einsteinstrasse 1 / 72172 Sulz a. N. / Germany

Tel.: +49 7454 9596-0 / Fax: +49 7454 9596-33 / e-mail: [sales@vbm-medical.de](mailto:sales@vbm-medical.de) / [www.vbm-medical.de](http://www.vbm-medical.de)

Connect with us:



CE / CE 0123

FLY016\_3.0\_EN